

Cannulas & Needles

Sterimedix Silkann Injection Cannulas




Europe's leading cannula manufacturer for the injection of fat and dermal fillers. Rigid and ultra-thin wall design. Made in Britain for 30 years to the highest possible quality standards. Gold standard used by leading clinical and industry opinion leaders. All of our cannulas are conveniently packed with a pre-hole needle.

Sterimedix Limited
T: +44 (0)1527 501480
E: sales@sterimedix.com
www.sterimedix.com

Sterimedix
aesthetic

Cannulas & Needles

Sterimedix Silkann GTI Cannula



Europe's leading cannula manufacturer bring you a unique GROOVED-TIP cannula for small scars and surface defects of the face designed by a renowned plastic surgeon in the field of fat transfer. Safely divide tethered fibrous bands of scar to restore contour defects.

Sterimedix Limited
T: +44 (0)1527 501480
E: sales@sterimedix.com
www.sterimedix.com

Sterimedix
aesthetic

Compression Garments

Triangle Surgical Formex Garments



Soft and comfortable compression garments with true four-way stretch material for optimum results. USA made.

- Sensibly priced, cost saving
- Huge range of sizes and styles in stock
- New styles design service available

Triangle Surgical Ltd
T: +44 (0)1494 876990
E: info@trianglesurgical.uk
www.trianglesurgical.uk

Dermatoscopes

Heine DELTA 30 Dermatoscope



Reliable diagnostics due to our best imaging
Heine's new high-end dermatoscope; standing out due to the largest true field of view.

- Fast change from contact to non-contact examination
- One button polarisation & non-polarised
- Brightness control
- 10x magnification
- USB rechargeable or optional desk charger
- Smartphone compatible

SurgiSol Ltd
T: 0845 508 4651
E: info@surgisol.com
surgisol.com

S
SURGISOL

Dermatoscopes

Heine DELTAone



Small on the outside, big on the inside
A Dermatoscope uniting diagnostic precision & mobility

– Crystal-clear images are accomplished using the new achromatic HEINE optical system.

- Switch between polarised & non-polarised light at the touch of a button
- Brightness control
- Easily cleaned
- 10x magnification
- USB rechargeable
- Smartphone compatible

SurgiSol Ltd
T: 0845 508 4651
E: info@surgisol.com
surgisol.com

S
SURGISOL

Energy Devices - Lasers

Aerolase Era



Aerolase Era is the only short-pulse Erbium:YAG laser. The 350 microsecond (350µs) pulse is shorter than the thermal relaxation time of skin, $\approx 800\mu s$, which means there is no thermal build up in the skin, minimal damage to surrounding tissue making recovery time days rather than weeks for full (not fractional) laser resurfacing.

Belle
T: +44 (0)1803 298415
E: mail@belle.org.uk
www.belle.org.uk

BELLE
FOR A LIFE-CHANGING RECOVERY